

# Microvessel Analysis


Detect new angiogenic vessels from larger vasculature, even in overstained xenografts.

Analyzing capillary and microvessel infrastructure is critical in many disease areas. From matrigel assays in cardiovascular research to vascular morphology changes in xenografts to microcapillary measurements in ophthalmology, researchers need to understand and quantify vascular architecture. With whole slide scanning and the Microvessel Analysis Algorithm, a tedious and subjective manual counting task can be automated and standardized.

Automate and standardize vascular architecture measurements with microvessel analysis across whole slides.

Microvessel analysis typically starts with a Factor VIII or CD31 stained slide, scanned at 20x or 40x. A novel dual-level thresholding algorithm screens out non-specific staining. Regions are then assembled into vessels, and individual statistics are calculated on each vessel. The user can choose which statistic to graph, and then can generate a summary report.

Statistic	Tumor	Normal
Number of Vessels	704	352
Total Analysis Area ( $\mu\text{m}^2$ )	12811307	4110861
Total Stain Area ( $\mu\text{m}^2$ )	9535109	62105
Average Stain Intensity	158667	191
Mean Vessel Area ( $\mu\text{m}^2$ )	143066	213
Median Vessel Area ( $\mu\text{m}^2$ )	996743	121
Standard Deviation of Vessel Area	13516	302
Mean Vessel Perimeter ( $\mu\text{m}$ )		
Median Vessel Perimeter ( $\mu\text{m}$ )		
Standard Deviation of Vessel Perimeter ( $\mu\text{m}$ )		
Mean Lumen Area ( $\mu\text{m}^2$ )		
Median Lumen Area ( $\mu\text{m}^2$ )		
Standard Deviation of Lumen Area ( $\mu\text{m}^2$ )		
Mean Vascular Area ( $\mu\text{m}^2$ )		
Median Vascular Area ( $\mu\text{m}^2$ )		
Standard Deviation of Vascular Area ( $\mu\text{m}^2$ )		
Mean Vessel Wall Thickness ( $\mu\text{m}$ )		
Median Vessel Wall Thickness ( $\mu\text{m}$ )		
Standard Deviation of Vessel Wall Thickness ( $\mu\text{m}$ )		



## Microvessel Analysis Algorithm

Microvessel Analysis is a smart object recognition algorithm that can filter out existing larger vessels, and identify and compute morphology statistics on only the newer angiogenic vessels.

# Specifications

## INPUT PARAMETERS

- Overlap Size
- Image Zoom
- Markup Image Type
- Include/Exclude Vessels Without Lumen
- Filtering/Smoothing Level
- Dark Staining Threshold
- Light Staining Threshold
- Region Joining Parameter ( $\mu\text{m}$ )
- Vessel Completion Parameter ( $\mu\text{m}$ )
- Minimum Vessel Area Threshold ( $\mu\text{m}^2$ )
- Maximum Vessel Area Threshold ( $\mu\text{m}^2$ )
- Maximum Vessel Wall Thickness ( $\mu\text{m}$ )
- Output Histogram Statistic
- Histogram Start Value ( $\mu\text{m}$  or  $\mu\text{m}^2$ )
- Histogram End Value ( $\mu\text{m}$  or  $\mu\text{m}^2$ )
- Number of Bins
- Endothelial Stain (RGB)
- Background Stain (RGB)
- 2nd Background Stain (RGB)

## SUMMARY OUTPUT STATISTICS

- Number of Vessels
- Total Analysis Area ( $\mu\text{m}^2$ )
- Total Stain Area ( $\mu\text{m}^2$ )
- Average Stain Intensity
- Microvessel Density (vessels/ $\mu\text{m}^2$ )
- Mean, Median, SD Vessel Area ( $\mu\text{m}^2$ )
- Mean, Median, SD Vessel Perimeter ( $\mu\text{m}$ )
- Mean, Median, SD Vascular Area ( $\mu\text{m}^2$ )
- Mean, Median, SD Vessel Wall Thickness ( $\mu\text{m}$ )
- Average Optical Density (RGB)

## INDIVIDUAL VESSEL STATISTICS

- Vessel Area ( $\mu\text{m}^2$ )
- Vessel Perimeter ( $\mu\text{m}$ )
- Lumen Area ( $\mu\text{m}^2$ )
- Vascular Area ( $\mu\text{m}^2$ )
- Vessel Wall Thickness ( $\mu\text{m}$ )

## MARK-UP TYPES

- 0 – Tuning – Vessel or Endothelial Stain
- 1 – Tuning – Vessel Thresholding
- 2 – Analysis – Vessels
- 3 – Analysis – Lumen and Vascular Cells
- 4 – Analysis – Vessels (Show Excluded Vessels)

**HEADQUARTERS**  
 Aperio  
 1360 Park Center Drive  
 Vista, CA 92081

**EUROPEAN OPERATION**  
 Aperio  
 3 The Sanctuary  
 Eden Office Park  
 Ham Green  
 Bristol BS20 0DD, UK

Main: 866.478.4111  
 Fax: 760.539.1116

Main: 011 44 1275 375 123  
 Fax: 011 44 1275 373 501

info@aperio.com  
 support@aperio.com  
 www.aperio.com

europa@aperio.com  
 europe@aperio.com  
 www.aperio.com

Aperio products are FDA cleared for specific clinical applications, and are intended for research use for other applications.

Copyright © 2008, Aperio. All rights reserved. ScanScope® and ImageScope™ are trademarks of Aperio.

IMC-1051-03132008

## MARK-UP TYPES




Endothelial cells




Endothelial cells grouped to make regions




Vessel perimeter drawn in middle of vessel walls



The algorithm allows the user to specify the distance at which two regions should be joined (RJ), and the distance for which a vessel should be completed (VC).



Calculate vessel wall thickness changes across whole slide images.



Microvessel Analysis can quantify vascular changes extending outward from a xenograft, and then compare with normal tissue vasculature.

For more information about our solutions for digital pathology, please visit [www.aperio.com](http://www.aperio.com) or call us at 866-478-4111.



bringing digital pathology to life